Silene™  
Peripheral Covered Stent System  
**INSTRUCTIONS FOR USE**

ENGLISH

Silene™ Peripheral Covered Stent System

**INSTRUCTIONS FOR USE**

Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.

**1.0 DEVICE NAME**

The device brand name is Silene™ Peripheral Covered Stent System.

**2.0 DESCRIPTION**

The Silene™ Peripheral Covered Stent System is pre-mounted onto a balloon expandable delivery catheter and has an ePTFE covering. The principal mechanism of action for the Silene™ Peripheral Covered Stent is inflation to the balloon via an inflation device.

The delivery catheter is compatible with 0.035" (0.89mm) guidewire and has a usable length of between 120-135cm. The delivery catheter has 2 radiopaque (RO) marker bands on either side of the crimped stent to aid in stent placement and positioning under fluoroscopy.

The Silene™ Peripheral Covered Stent System is packaged sterile and intended for single use only.

Stent Platform:

The stent is provided in a range of sizes from diameters 6.0–12.0mm, and lengths of 15-70mm. The Silene™ is based on the United States Patented intravascular stent design – DIRECT-STENT® intravascular Stent design. The alloy is medical grade L605 Cobalt Chromium Alloy.

Delivery Catheter:

The Over the Wire Delivery (OTW) Catheter has a balloon located at its distal end designed for peripheral transluminal angioplasty (PTA) by percutaneous way. The catheter is a dual-lumen shaft. One lumen is used for inflation of the balloon and the other permits the use of a guidewire (0.035") to facilitate advancement of the stent system to and through the stenosis to be dilated. The balloon has two radiopaque markers to aid positioning the balloon relative to the stenosis. The material used for the balloon allows the operator to control the external diameter depending on the inflation pressure. Refer to table with compliance chart for ATM pressure versus Balloon diameter. A rigid connector hub is used at the end of the proximal part. This hub has a female luer lock connector enabling connection with other medical devices to permit the inflation of the balloon with an inflation device.

Polymer:

The ePTFE polymer consists of 100% virgin PTFE homopolymer. The polymer is manufactured from resin into a membrane that is wrapped to form a tube. The method of manufacturing produces an inherently high strength ePTFE covering with thin wall. The porosity of the polymer and the lattice structure is proprietary including the manufacturing process.

**SIZE MODELS**

PRODUCT DESCRIPTION (mm)	REFERENCE NUMBER
Silene™ 6.0 x 15	756015
Silene™ 6.0 x 25	756025
Silene™ 6.0 x 35	756035
Silene™ 6.0 x 50	756050
Silene™ 6.0 x 70	756070
Silene™ 7.0 x 15	757015
Silene™ 7.0 x 25	757025
Silene™ 7.0 x 35	757035
Silene™ 7.0 x 50	757050
Silene™ 7.0 x 70	757070
Silene™ 8.0 x 15	758015
Silene™ 8.0 x 24	758024
Silene™ 8.0 x 32	758032
Silene™ 8.0 x 50	758050
Silene™ 8.0 x 70	758070
Silene™ 9.0 x 15	759015
Silene™ 9.0 x 24	759024
Silene™ 9.0 x 32	759032
Silene™ 9.0 x 50	759050
Silene™ 9.0 x 70	759070
Silene™ 10.0 x 15	7510015
Silene™ 10.0 x 23	7510023
Silene™ 10.0 x 31	7510031
Silene™ 10.0 x 48	7510048
Silene™ 10.0 x 65	7510065
Silene™ 12.0 x 14	7512014
Silene™ 12.0 x 23	7512023
Silene™ 12.0 x 30	7512030
Silene™ 12.0 x 43	7512043
Silene™ 12.0 x 58	7512058

**3.0 CONTENTS**

- 1 Silene™ Peripheral Covered Stent System
- 1 Instructions for Use
- 1 Patient Implant Card

**SUPPLIED:**

- STERILE
- Sterilized with Ethylene Oxide gas (EtO)
- FOR SINGLE USE ONLY
- DO NOT RESTERILIZE

**STORAGE:**

- Keep from Sunlight
- Keep Dry
- Do not use if the package has been opened or damaged
- Use by the "Use By" date noted on the package
- Use only between 15 °C and 30 °C

**OTHER:**

- Indicates the date when the medical device was manufactured
- Indicates the medical device manufacturer
- Indicates the manufacturer's batch code for identification
- Indicates the manufacturer's catalogue number for identification
- Items marked MR conditional pose no known hazards in specified MR environments with specified conditions for use
- Quantity

#### 4.0 INDICATIONS FOR USE / INTENDED USE:

The SILENE™ Peripheral Covered Stent System is indicated for restoring and improving the patency of the iliac artery.

#### 5.0 CONTRAINDICATIONS:

The Silene™ Peripheral Covered Stent is contraindicated for use in:

- 5.1 Tortuous anatomy with difficulty accessing the site with standard stent placement techniques.
- 5.2 With a lesion judged to prevent complete inflation of an angioplasty balloon.
- 5.3 Aneurysms adjacent to the site of stent implantation.
- 5.4 In whom antiplatelet and/or anticoagulation therapy is contraindicated.
- 5.5 Patients who exhibit angiographic evidence of existing fresh thrombus.
- 5.6 Patients with a contraindication for anti-platelet/anti-coagulant therapy. This includes patients who have had major surgery, an obstetrical delivery, organ biopsy, or puncture of a non-compressible vessel within 14 days of this procedure. Also excluded are those patients with a history of gastrointestinal bleeding, recent C.V.A. diabetic hemorrhagic retinopathy, or any other condition compromised by prolonged anti-coagulation.
- 5.7 Target lesions in an area involving a major side branch. A major side branch is defined as a vessel that would be bypassed if it were occluded.
- 5.8 Patients who have experienced a recent (less than 1 week) acute myocardial infarction.
- 5.9 To reduce the potential for vessel damage, the inflated diameter of the balloon should approximate the diameter of the vessel just proximal and distal to the stenosis.
- 5.10 Stenting in patients who are not acceptable candidates for peripheral artery bypass surgery requires careful consideration, including possible hemodynamic support during stenting, as treatment of this patient population carries special risk.
- 5.11 The use of adjacent stents, or in-stenting, or overlapping of stent, or anchoring to another stent is contraindicated.

#### 6.0 POTENTIAL ADVERSE EFFECTS/COMPLICATIONS

- 6.1 Death
- 6.2 Arrhythmia
- 6.3 Peripheral Artery Bypass Surgery
- 6.4 Abscess
- 6.5 Drug reactions to antiplatelet agents/contrast medium
- 6.6 Pyrogenic reactions
- 6.7 Hypotension
- 6.8 Inadequate stent graft implantation or intimal trauma
- 6.9 Intimal flap
- 6.10 Stent misplacement, migration, or deformation
- 6.11 Endocarditis
- 6.12 Hemorrhage, requiring transfusion
- 6.13 Hypotension/Hypertension
- 6.14 Infection and pain at insertion site
- 6.15 Ischemia, myocardial
- 6.16 Perforation, dissection, injury, or rupture of the vessel.
- 6.17 Pseudoaneurysm, femoral
- 6.18 Restenosis of stented segment
- 6.19 Vessel Spasm
- 6.20 Vessel thrombosis
- 6.21 Stroke/Cerebrovascular Accidents
- 6.22 Total occlusion of peripheral artery
- 6.23 Sepsis/infection

#### 7.0 WARNINGS AND CAUTIONS

The implantation of the stent should be performed only by physicians who have received adequate training in interventional techniques in stent placement. The following instruction provide technical guidance, but do not obviate the need for formal training in the use of this device.

- 7.1 Carefully inspect the stent system prior to use and ensure it has not been damaged during shipment.
  - 7.2 Ensure that the stent dimensions is suitable for the specific procedure.
  - 7.3 Do not manipulate the stent covering, remove, reposition, or hand crimp the covered stent.
  - 7.4 Do not use if the inner pouch is opened or damaged.
  - 7.5 Do not force advancement or withdrawal of the guidewire or catheter delivery system if resistance is encountered.
  - 7.6 Do not use power injection system to inflate the balloon catheter. Use only hand-held inflation device specifically design for stent procedures.
  - 7.7 Do not use air, or other gaseous medium to inflate the balloon.
  - 7.8 Patients with allergies to Cobalt Chromium L605 alloy or PTFE may suffer an allergic response to this implant.
  - 7.9 Prior to expanding the stent, utilize fluoroscopy to verify that the stent has been properly positioned. Do not deploy the stent if it has not been properly positioned.
  - 7.10 For single use only. Do not reuse, reprocess, or resterilize as this may compromise the stent covering, structural integrity, and result in patient injury, illness, or death.
  - 7.11 Balloon pressures are indicated on the product label and should be monitored during inflation. To ensure full stent expansion, inflate to at least the nominal pressure. Do not exceed Rated Burst Pressure.
- WARNING: For the maximum recommended post dilation balloon pressure and diameter see the table below. Post dilation greater than the recommended diameter may result in the structural damage of the stent causing serious injury or death. Under-expansion of the stent may result in stent migration.

system is unable to reach or cross the lesion easily, the procedure should be aborted. Refer to the section "Removal of an unexpanded stent".

- 7.12 Do not attempt to attempt to reposition a fully or partially expanded stent.
- 7.13 Do not pull an unexpanded stent back through a sheath or a guiding sheath.
- 7.14 Do not stent across in a bifurcated vessel.
- 7.15 Do not overlap stent manufactured from different materials.
- 7.16 Do not implant device in patients who have an active infection at the implant site.

#### 8.0 Preparation:

For preparation of the Silene™ Peripheral Covered Stent System, it is recommended that in addition to the stent delivery catheter, the following equipment be assembled for use.

- 8.1 Guidewire (0.035")
- 8.2 Guidewire insertion tool -Optional
- 8.3 Balloon inflation device (i.e. an inflation device with a stopcock if desired and connector tubing)
- 8.4 Sterile bowl with heparinized flush solution
- 8.5 Sterile bowl with a 1:1 mixture of low viscosity contrast medium and sterile saline
- 8.6 10cc or larger plastic luer lock syringe (for stent delivery preparation)

#### 9.0 GENERAL PROCEDURE

- 9.1 Patient preparation should be the same as for any angioplasty procedure.
- 9.2 Standard techniques for placement of a sheath and guidewire should be used when using this device.

CAUTION: Sheath size should be used as labeled on the box and inside pouch

- 9.3 Perform diagnostic angiography to confirm site of implantation and branch vessels.
- 9.4 Evaluate and mark the site of implant.

CAUTION: Standard PTA techniques should be used if predilation is necessary.

#### 10.0 STENT SIZING

- 10.1 Measure the length of the lesion to determine the stent required. The Silene™ Peripheral Covered Stent is labeled to expand to the diameter and length. In some cases, due to anatomical challenges it may foreshorten slightly more than as labeled. It is recommended to have stent extended at least 5mm proximal and distal to the lesion after expansion of the stent when possible.
- 10.2 Measure the inner diameter of the vessel to determine the approximate diameter stent to implant. The Silene™ is labeled as the outer diameter.

#### 11.0 Instructions For Use

CAUTION: Exercise care during device handling to reduce the possibility of disrupting the delicate placement of the stent on the balloon and accidental breakage, bending, or kinking of the stent delivery system shaft.

- 11.1 Gently remove the tray from the protective pouch.
- 11.2 Gently remove the stent system from the tray.

CAUTION: Forceful movements may dislocate the protector and the stent.

- 11.3 Carefully remove the balloon/stent protector sheath

CAUTION: Always pull at the very distal end of the protector sheath to avoid dislocation of the stent.

- 11.4 Visually check the stent crimping for uniformity, no protruding struts, and centering on the balloon between the 2 marker bands.
- 11.5 Connect a syringe containing sterile saline to the guidewire lumen and carefully flush the guidewire lumen.
- 11.6 Remove the syringe.

WARNING: Use only an appropriate balloon inflation medium (e.g. 50:50 mixture by volume of contrast medium and saline). Never use air or any gaseous medium to inflate the balloon.

- 11.7 Attached the pre-filled inflation device to the catheter hub for balloon inflation.

CAUTION: Ensure no bubbles are in the inflation device.

- 11.8 Ensure that the stent is in between the RO markerbands to verify all preparation steps did not create positive pressure which could cause inadequate stent deployment and the stent remains securely positioned within the RO markers prior to stent implant.
- 11.9 Leave the prepared stent system at ambient pressure.

CAUTION: DO NOT apply negative pressure to the delivery system prior to placement of the stent across the lesion. This may cause loosening and dislodgement of the stent.

#### 12.0 DEPLOYMENT

- 12.1 Position the guidewire under fluoroscopy in accordance with PCI techniques for stent implant.
  - 12.2 Back-load the proximal end of the guidewire into the distal tip of the delivery system until it exits at the guidewire exit port.
  - 12.3 Carefully insert the stent system through the Sheath over the guidewire.
- CAUTION: Sheath size should be used as indicated on the label.
- 12.4 Carefully advance the stent system into the peripheral artery over the guidewire while maintaining stable guidewire placement across the target lesion and under direct fluoroscopy visualization.
  - 12.5 Utilize fluoroscopy to verify the stent position within the lesion using the balloon radiopaque markers as reference points and assure adequate coverage of the lesion including proximal and distal margins.

CAUTION: DO NOT apply excessive force whilst accessing or crossing the lesion. This may damage the stent and/or dislodge it from the balloon. If resistance is felt at any time, stop the procedure, and determine the cause of resistance before proceeding. If the stent

- 12.6 Inflate the balloon gradually to expand the stent to the nominal deployment pressure as device label. Hold that pressure for 15-30 seconds.

CAUTION: Inflate to at least the nominal pressure (NP) indicated on the product label. DO NOT exceed Rated Burst Pressure (RBP).

- 12.7 Deflate the balloon by pulling vacuum on the inflation device to its maximum volume for more than 40 seconds in an upright position with the handle in up position.
- 12.8 If necessary, the delivery system balloon may be inflated once more to achieve optimum seating of the implanted stent.

CAUTION: Ensure no bubbles are in the inflation device.

- 12.9 If the stent is still not completely opposed to the vessel wall, the stent can be re-crossed and further expanded with a larger diameter balloon catheter.

NOTE: It is recommended that the guidewire remain across the lesion and stent implant site until the procedure is completed.

CAUTION: DO NOT post-dilate the stent to more than the maximum expandable diameter recommended in the compliance chart.

CAUTION: Exercise care when crossing a newly deployed stent with a peripheral guidewire, balloon catheter or stent system to avoid disruption of the stent geometry and stent migration.

**13.0 BALLOON DEFLATION AND DELIVERY SYSTEM REMOVAL**

- 13.1 Deflate the balloon in accordance with standard PCI procedures. Apply negative pressure to the balloon for at least 40 seconds. Before carefully pulling the delivery system out of the vessel.
- 13.2 If the balloon cannot be withdrawn from the stent easily, slightly advance and retract the delivery system very carefully until it is possible to pull it out.
- 13.3 Under direct fluoroscopic visualization retrieve the delivery system carefully into the sheath or guiding sheath.

CAUTION: If resistance is encountered while doing so, remove the delivery system and the guiding catheter as a single unit according to the instructions below.

**14.0 REMOVAL OF AN EXPANDED STENT**

- 14.1 Extreme care must be used when removal of an unexpanded or partial expanded stent is necessary.
- 14.2 Slowly pull back the stent system. The proximal end of the stent must be aligned with the distal tip of the sheath. The sheath and the stent system should then be removed as one unit.

CAUTION: If any resistance is felt while doing so, the stent system and the guiding catheter must be removed as a single unit.

CAUTION: DO NOT RE-INSERT the Silene™ Covered Stent. The stent may dislodge.

WARNING: Do Not attempt to pull back the stent into the sheath.

**15.0 MAGNETIC RESONANCE IMAGING**

Non-clinical testing has demonstrated that the Covered stent is single and two overlapped versions) is MR Conditional. A patient with this device can be scanned safely in an MR system immediately after placement under the following conditions:

- 15.1 Static magnetic field of 1.5-Tesla and 3-Tesla, only.
- 15.2 Maximum spatial gradient magnetic field of 4,000-gauss/cm (40-T/m)(extrapolated) or less.
- 15.3 Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 2-W/kg for 15 minutes of scanning (i.e., per pulse sequence) in the Normal Operating Mode of operation for the MR system.
- 15.4 Under the scan conditions defined about the Silene Covered Stent is expected to produce a maximum temperature rise of 2.3°C after 15-minutes of continuous scanning (i.e., per pulse sequence).

**16.0 REFERENCES**

The physician should consult recent literature on current medical practice on peripheral stenting, such as that published by ACC/AHA.

**BALLOON COMPLIANCE CHART**

(ATM)	6.0mm	7.0mm	8.0mm	9.0mm	10.0mm	12.0mm	
<b>6</b>	<b>5.96</b>	<b>6.96</b>	<b>8.00</b>	<b>9.00</b>	<b>10.00</b>	-	<b>Nominal</b>
7	6.02	7.05	8.08	9.04	10.13	<b>12.11</b>	
8	6.05	7.09	8.08	9.14	10.28	12.38	
9	6.10	7.19	8.17	9.29	10.27	12.57	
10	6.22	7.28	8.30	9.42	10.46	<b>12.62</b>	<b>Rated Burst</b>
11	6.25	7.35	8.44	9.53	10.61		
12	6.37	7.45	8.58	9.63	10.86		
13	6.42	7.50	8.72	9.82	10.90		
<b>14</b>	<b>6.42</b>	<b>7.54</b>	<b>8.80</b>	<b>9.99</b>	<b>11.00</b>	<b>Rated Burst</b>	

Assure full deployment of the stent based in balloon compliance. Deployment pressures should be based on lesion characteristics.

InSitu Technologies, Inc.  
 539 Phalen Boulevard  
 St. Paul, MN 55130  
 United States of America  
 Tel 651 389 1017  
 Email: QA-RA@insitu-tech.com

FDA Registered Facility  
 ISO 13485 Certified  
 For Export Only

**DISCLAIMER**

THERE IS NO EXPRESS OR IMPLIED WARRANTY, INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE, ON THIS ITEM. PRODUCT(S) DESCRIBED IN THIS PUBLICATION UNDER NO CIRCUMSTANCES SHALL BE LIABLE FOR ANY DIRECT, INCIDENTAL OR CONSEQUENTIAL DAMAGES OTHER THAN AS EXPRESSLY PROVIDED BY SPECIFIC LAW. NO PERSON HAS THE AUTHORITY TO BIND INSITU TO ANY REPRESENTATION OR WARRANTY EXCEPT AS SPECIFICALLY SET FORTH HEREIN.

Descriptions or specifications in InSitu Technologies printed matter, including this publication, are meant solely to generally describe the product at the time of manufacture and do not constitute any express warranties.

# ITALIAN

## Silene™ Peripheral Covered Stent System

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### SUPPLIED:

STERILE

Sterilized with Ethylene Oxide gas (EtO)

FOR SINGLE USE ONLY

DO NOT RESTERILIZE

### STORAGE:

Keep from Sunlight

Keep Dry

Do not use if the package has been opened or damaged

Use by the "Use By" date noted on the package

Use only between 15°C and 30°C

### OTHER:

Indicates the date when the medical device was manufactured

Indicates the medical device manufacturer

Indicates the manufacturer's batch code for identification

Indicates the manufacturer's catalogue number for identification

Items marked MR conditional pose no known hazards in specified MR environments with specified conditions for use

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