

ASAHI
ULTIMATEebros **3**
PTCA GUIDE WIRE

ASAHI ULTIMATEbros 3


PTCA GUIDE WIRE

Первый проводник для реканализации хронической окклюзии

- Гидрофильное покрытие длиной 40 см обеспечивает высокую проходимость в сильно измененных участках артерии.
- Специальная конструкция кончика упрощает проведение в истинный просвет и сокращает риск расширения ложного просвета.

Информация для заказа

■ Структура



Product	Catalog No.	Diameter	Usable Length	Spring Coil Length	Radiopaque Length	Tip Shape
ULTIMATEbros 3	AHW14S003S	0.36mm (0.014inch)	180cm	11cm	11cm	Straight
	AHW14S303S	0.36mm (0.014inch)	300cm	11cm	11cm	Straight
	AHW14S003J	0.36mm (0.014inch)	180cm	11cm	11cm	J-Tip
	AHW14S303J	0.36mm (0.014inch)	300cm	11cm	11cm	J-Tip

SLIP COAT® is a trademark of Angiotech Biocoatings Corporation.

Manufactured by:



ASAHI INTECC CO.,LTD.

3-100 Akatsuki-cho, Seto, Aichi 489-0071 Japan
TEL: +81-561-48-5551 FAX: +81-561-48-5552
<http://www.asahi-intecc.com>

Distributed by:



www.raymed.ru

ASAHI INTECC CO., LTD. EUROPE

Strawinskylaan 967
(WTC Tower D-9th floor)
1077 XX Amsterdam,
The Netherlands
TEL: +31-(0)20-794-0640
FAX: +31-(0)20-794-0641

ASAHI INTECC USA, INC.

2500 Red Hill Avenue, Suite 210
Santa Ana, California,
92705, USA
TEL: +1-949-756-8901
FAX: +1-949-756-8165

ASAHI INTECC (THAILAND) CO., LTD.

158/1 Moo5,
Bangkadi Industrial Park,
Tiwanon Road, Tambol Bangkadi,
Amphur Muang,
Pathumthani 12000,
Thailand
TEL: +66-2-501-1302
FAX: +66-2-501-1305

ASAHI INTECC CO., LTD.
SINGAPORE BRANCH

3 Bishan Place
#06-01 Bishan CPF Building
Singapore 579838
TEL: +65-6250-4248
FAX: +65-6250-4247

ASAHI INTECC CO., LTD.
BEIJING REPRESENTATIVE OFFICE

Room 18, 10th Floor,
Twin Towers (East),
B12 Jianguomenwai Ave.,
Beijing, China 100022
TEL: +86-10-5123-5100
FAX: +86-10-5123-5177

ASAHI INTECC CO., LTD.
HONG KONG BRANCH

Unit 501-2, 5/F, Hon Kwok Jordan Centre
No 7 Hillwood Road,
Tsim Sha Tsui, Kowloon
Hong Kong
TEL: +852-2316-2510
FAX: +852-2316-2485