

ASAHI  
RG3

PTCA GUIDE WIRE

# ASAHI RG3

PTCA GUIDE WIRE


## Бережность к артериям. Легкая управляемость.

Безопасный и эффективный проводник для осуществления ретроградного доступа при реканализации коронарных артерий.

Оптимальная жесткость проводника, гидрофильное покрытие и 0.26 мм (0.010") профиль обеспечивают более высокую проходимость. Благодаря уменьшению риска повреждения внутренней стенки во время прохождения через извилистые участки сосудов, риск возникновения осложнений сведен к минимуму.

### Информация для заказа

#### ■ Структура



Product	Catalog No.	Diameter	Usable Length	Spring Coil Length	Radiopaque Length	Tip Shape
<b>ASAHI RG3</b>	AHW10S302S	0.26mm (0.010inch)	330cm	8cm	3cm	Straight

SLIP COAT® is a trademark of Angiotech Biocoatings Corporation.

Manufactured by:



**ASAHI INTECC CO.,LTD.**

3-100 Akatsuki-cho, Seto, Aichi 489-0071 Japan

TEL: +81-561-48-5551 FAX: +81-561-48-5552

<http://www.asahi-intecc.com>

Дистрибьютор:



**RAYMED CENTER**

[www.raymed.ru](http://www.raymed.ru)

ASAHI INTECC CO., LTD. EUROPE

Strawinskylaan 967  
(WTC Tower D-9th floor)  
1077 XX Amsterdam,  
The Netherlands  
TEL: +31-(0)20-794-0640  
FAX: +31-(0)20-794-0641

ASAHI INTECC USA, INC.

2500 Red Hill Avenue, Suite 210  
Santa Ana, California,  
92705, USA  
TEL: +1-949-756-8901  
FAX: +1-949-756-8165

ASAHI INTECC (THAILAND) CO., LTD.

158/1 Moo5,  
Bangkadi Industrial Park,  
Tiwanon Road, Tambol Bangkadi,  
Amphur Muang,  
Pathumthani 12000,  
Thailand  
TEL: +66-2-501-1302  
FAX: +66-2-501-1305

ASAHI INTECC CO., LTD.  
SINGAPORE BRANCH

3 Bishan Place  
#06-01 Bishan CPF Building  
Singapore 579838  
TEL: +65-6250-4248  
FAX: +65-6250-4247

ASAHI INTECC CO., LTD.  
BEIJING REPRESENTATIVE OFFICE

Room 18, 10th Floor,  
Twin Towers (East),  
B12 Jianguomenwai Ave.,  
Beijing, China 100022  
TEL: +86-10-5123-5100  
FAX: +86-10-5123-5177

ASAHI INTECC CO., LTD.  
HONG KONG BRANCH

Unit 501-2, 5/F, Hon Kwok Jordan Centre  
No 7 Hillwood Road,  
Tsim Sha Tsui, Kowloon  
Hong Kong  
TEL: +852-2316-2510  
FAX: +852-2316-2485