


ANGIO-LINE™

Products and Accessories for Angiography



With extensive experience and proven track record in angiography and cardiology, Biometrix offers a complete product line for interventional cardiology and radiology. Contrast media injection lines, sets and syringes, angiography guide wires, haemostatic valve introducer sets, PTCA accessories kits and comprehensive custom procedure trays – part of *BIOMETRIX CPT* program – are some of the Angio-Line™ products and accessories.



Hemostatic Valve Introducers

Standard Length Introducers

Code	Dilator Size	Needle	Packaging	Qty./Box
HD-6301	5Fr	+	Tray	10
HD-6401	6Fr	+	Tray	10
HD-6501	7Fr	+	Tray	10
HD-6701	8Fr	+	Tray	10
HD-6801	9Fr	+	Tray	10
HD-7301	5Fr	-	Pouch	20
HD-7401	6Fr	-	Pouch	20
HD-7501	7Fr	-	Pouch	20
HD-7701	8Fr	-	Pouch	20
HD-7901	9Fr	-	Pouch	20
HD-X6301	5Fr	+	Pouch	50
HD-X6401	6Fr	+	Pouch	50
HD-X6501	7Fr	+	Pouch	50
HD-X6701	8Fr	+	Pouch	50
HD-X6801	9Fr	+	Pouch	50

Set Contents - Standard Length

- Sheath Dilator Hemostatic Valve Unit with Side Arm and Stopcock
- 0.038" x 50cm Double Ended Guidewire
- 18Ga x 7cm Introducer Needle (in needle sets)

Short Length Introducers

Code	Dilator Size	Needle	Packaging	Qty./Box
HD-P6201	4Fr	+	Tray	10
HD-P6301	5Fr	+	Tray	10

Micropuncture Introducer Sets

Standard Length Introducers

Code	Dilator Size
HD-Y6201	4Fr
HD-Y6301	5Fr
HD-Y6401	6Fr

Qty/Box: 10

Set Contents

- 13 cm Sheath Dilator Hemostatic Valve Unit with Side Arm and Stopcock
- 0.018" x 40cm Stainless Steel Guidewire
- 21Ga x 3.8cm Introducer Needle with Echogenic Tip
- Syringe 10ml

Long Length Introducers

Code	Dilator Size
HD-X9301	5Fr
HD-X9401	6Fr
HD-X9501	7Fr

Qty/Box: 10

Set Contents

- 23 cm Sheath Dilator Hemostatic Valve Unit with Side Arm and Stopcock
- 0.018" x 70cm Stainless Steel Guidewire
- 21Ga x 7cm Introducer Needle with Echogenic Tip



Long Length Introducers

Code	Dilator Size
HD-9301	5Fr
HD-9401	6Fr
HD-9501	7Fr
HD-9701	8Fr

Qty/Box: 10

Set Contents - Long Length

- 23 cm Sheath Dilator Hemostatic Valve Unit with Side Arm and Stopcock
- 0.038" x 70cm Double Ended Guidewire
- 18Ga x 7cm Introducer Needle

Set Contents - Short Length

- 7 cm Sheath Dilator Hemostatic Valve Unit with Side Arm and Stopcock
- 0.018" x 40cm Double Ended Guidewire
- 21Ga x 3.8cm Introducer Needle

Power Syringes

Code	Description	Qty/Box
BI-0538	150ml Injection Power Syringe with Filling Tube for Medrad Mark V and Mark V Ius Injectors	50
BI-0540	150ml Power Syringe for Angiomat 6000 Injector with Filling Tube	60
BI-0548	200ml Injection Power Syringe with Filling Tube for Medrad Mark V and Mark V Plus Injectors	50
BI-0910	Filling Tube for Power Syringe	150



High Pressure Lines (I200psi)

Code	Dimension (mm)	Length	Material
BL-1075	1.9 x 4.7	75cm	PVC
BL-1125	1.9 x 4.7	125cm	PVC
BL-1160	1.9 x 4.7	160cm	PVC
BM-4050	1.8 x 3.7	50cm	Braided
BM-4075	1.8 x 3.7	75cm	Braided
BM-4120	1.8 x 3.7	120cm	Braided
BM-4150	1.8 x 3.7	150cm	Braided
BM-4160	1.8 x 3.7	160cm	Braided

Qty/Box: 50



- All lines feature rotating adaptor

Pressure Monitoring Lines

Code	Description	Qty/Box
AH-1015	15cm M/F	150
AH-1030	30cm M/F	150
AH-1050	50cm M/F	150
AH-1100	100cm M/F	150
AH-1150	150cm M/F	100
AH-1180	180cm M/F	100
AH-1200	200cm M/F	100
AH-2100	100cm M/F	150
AH-2150	150cm M/F	100
AH-2200	200cm M/F	100
AH-2240	240cm M/F	100

Qty/Box: 10



Guide Wires

Code	Diameter	Length	J Tip	Features
HR-3031	0.035"	180cm	3mm	FC, FS, DD
HR-3812	0.035"	150cm	-	FC, HB, STR
HR-3822	0.035"	150cm	3mm	FC, FS, HB
HR-3823	0.035"	150cm	3mm	FC, FS, FB
HR-3828	0.035"	150cm	3mm	FC, FS, HB, FB
HR-3831	0.035"	150cm	3mm	FC, FS, DD
HR-3922	0.035"	260cm	3mm	FC, FS, HB
HR-4822	0.038"	150cm	3mm	FC, FS, HB
HR-4862	0.038"	150cm	-	FC, HB, HD, STR
HR-4912	0.038"	260cm	-	FC, HB, STR
HR-4922	0.038"	260cm	3mm	FC, FS, HB
HR-5521	0.035"	150cm	3mm	MC

Qty/Box: 25

- FC - Fixed Core
- FS - Finger Straightenable
- HB - Heparin Bonded
- FB - Firmed Body
- HD - Heavy Duty
- DD - Double Distal
- STR - Straight
- MC - Movable Core



Drainage Catheter Connectors

Code	Description	Stopcock
DQ-0305	5cm	-
DQ-0325	25cm	-
DQ-0430	40cm	+

Qty/Box: 100



Angiography Waste Bag

Code	Description
BQ-0215	Single Line Waste Bag
BQ-0205	Double Line Waste Bag

Qty/Box: 25



Split Sheath Introducer Sets

Code	Dilator Size
JF-0701	7Fr
JF-0801	8Fr
JF-0901	9Fr
JF-0951	9.5Fr
JF-1001	10Fr
JF-1051	10.5Fr
JF-1101	11Fr

Qty/Box: 10

Set Contents

- 15cm Split Sheath/Dilator Assembly
- 0.038" x 60cm Double Ended Guidewire
- 18Ga x 7cm Introducer Needle
- 12cc Syringe LL



Micropuncture Sets

Code	Description	Guide Wire
HT-0210	4Fr Set	Nitinol
HT-0220	4Fr Set	Stainless Steel
HT-0310	5Fr Set	Nitinol
HT-0320	5Fr Set	Stainless Steel

Qty/Box: 10

Set Contents

- 4 Fr or 5 Fr Vessel Dilator Set
- 0.018" x 40cm Guidewire
- 21Ga x 7cm Introducer Needle with Echogenic Tip



Micropuncture Guidewire

Code	Diameter	Length	Tip	J Radius	Features
HR-2242	0.018"	40cm	Platinum Floppy 3 cm	3mm	Nitinol SS

Qty/Box: 20

Introducer Needles

Code	Cannula Size	No. of parts	Bevel	Wing
HQ-1510	18Ga x 7cm	One Part	Standard	-
HQ-1511	18Ga x 7cm	One Part	Standard	+
HQ-1514	18Ga x 7cm	One Part	Short	-
HQ-1521	18Ga x 7cm	Two Part	Standard	-
HQ-1525	18Ga x 9cm	One Part	Standard	-
HQ-1611	18Ga x 7cm	One Part	Standard	+

Qty/Box: 100



Finger Grip Syringe

Code	Description
BI-0121	12cc with Rotator

Qty/Box: 50



PTCA Accessories Kits

Code	Dilator Size
BQ-0805	PTCA Kit with Y-Angioplasty Valve
BQ-0905	PTCA Kit with Click Angioplasty Valve
BQ-1905	PTCA Kit with Click Angioplasty Valve and Extension Line
BQ-9050	PTCA Kit with Inflation Device

Qty/Box: 25, BQ-9050: 12



Set Contents for BQ-0805

- Y-Angioplasty Valve
- Torque Device
- Insertion Tool

Set Contents for BQ-0905

- Click Angioplasty Valve
- Torque Device
- Insertion Tool

Set Contents for BQ-1905

- Click Angioplasty Valve with 20 cm Side Arm and Stopcock
- Torque Device
- Insertion Tool

Set Contents for BQ-9050

- Inflation Device
- Click Angioplasty Valve
- Torque Device
- Insertion Tool

Inflation Device

Code	Description
BQ-9001	Inflation Device 20cc

Qty/Box: 10



Angioplasty Valves

Code	Description
BQ-0504	Single Y-Angioplasty Valve
BQ-0505	Click Angioplasty Valve
BQ-0602	Single Y-Angioplasty Valve with Insertion Tool
BQ-0701	Double Y-Angioplasty Valve
BQ-0702	Double Y-Angioplasty Valve with Insertion Tool

Qty/Box: 100



Torque Device

Code	Description
BQ-0010	Torque Device Guidewires 0.014"-0.022"

Qty/Box: 100



Angiography Manifolds and Stopcocks

High Pressure Manifolds			
Code	Ports & Position	Pressure (PSI)	Tube
BE-2110	2 R OFF	600	-
BE-2111	2 R OFF	600	15cm High Pressure
BE-2120	2 L OFF	600	-
BE-3110	3 R OFF	600	-
BE-3111	3 R OFF	600	15cm High Pressure
BE-3120	3 L OFF	600	-
BE-3210	3 R OFF	600	-



High Pressure Stopcock	
Code	Ports & Position
BE-1111	3 way High Pressure Stopcock with 10cm Pressure Line
BE-1149	3 way High Pressure Stopcock with Rotating Adaptor MLL

Qty/Box: 50

Contrast Media Manifold Sets

Code	Description
BQ-1000	Angiography Manifold 2 Port Set

Qty/Box: 25

Set Contents

- Infusion Set Vented
- Manifold 2 Port F/F with 15cm High Pressure Line with Rotating Adaptor
- 12ml Finger Grip Syringe



Code	Description
BQ-3000	Angiography Manifold 3 Port Set

Qty/Box: 25

Set Contents

- Infusion Set
- 3 Port Manifold F/Rot. MLL
- 100cm Infusion Line F/MLL
- 12ml Finger Grip Syringe

Code	Description
BQ-2000	Contrast Media Set

Qty/Box: 25

Set Contents

- Pressure Line
- Double Check Valve
- 3 Way High Pressure Stopcock
- Fluid Administration Set Vented
- 12ml Finger Grip Syringe

Coiled Low Pressure Sets

Code	Description
BO-1150	Coiled Low Pressure Connecting Tube M/F
BO-8150	Coiled Low Pressure Connecting Tube with One-Way Check Valve
BO-9150	Coiled Low Pressure Connecting Tube with T-Type Check Valve

Qty/Box: 25



Image and Protective Covers

Code	Description	Cover
BY-1001	R 65 cm	Image
BY-1002	80 x 80 cm	Image
BY-1003	100 x 100 cm	Image
BY-2001	80 x 100 cm	Protective
BY-2002	100 x 100 cm	Protective
BY-2003	80 x 110 cm	Protective
BY-2004	120 x 120 cm	Protective

Qty/Box: 50



Reinforced Angiography Drape

Code	Description
BY-3001	220 x 300cm One Sided Polyethylene Clear Film
BY-3002	240 x 300cm Two Sided Polyethylene Clear Film
BY-3003	240 x 400cm Two Sided Polyethylene Clear Film
BY-3010	240 x 300cm Heart Pacemaker Drape
BY-3011	240 x 225cm Heart Pacemaker Drape

Qty/Box: 20



Surgical Gowns

Code	Specifications	Material	Size
DY-2001	Reinforced	SMS	M
DY-2002	Reinforced	SMS	L
DY-2003	Reinforced	SMS	XL
DY-2004	Standard	SMS	M
DY-2005	Standard	SMS	L
DY-2006	Standard	SMS	XL
DY-2007	Standard	SMS	L Long
DY-2008	Standard	SMS	XL Long

Qty/Box: 25



Enhance your cath lab efficiency
using *BIOMETRIX CPT*

- Fully customized packs
- Saves time and improves staff efficiency
- Components from one single source
- Simplifies logistic chain
- Minimizes stock holding



Paardeweide 36 P.O.B. 7191 | 4800 GD Breda, The Netherlands
Tel: 31-76-5422720 | Fax: 31-76-5422716
sales@biometrixmedical.com | www.biometrixmedical.com