




Angiographic Kits



Bioptimal™

Angiographic Kits

Product Description

Bioptimal Angiographic Kits consists of one 3-gang manifold, one IV administration set, one contrast media administration set, one pressure tubing and one control syringe.

The male rotating connector at the distal end of the manifold will be connected to PTCA catheters or angiographic catheters that will be inserted into the patient using Seldinger techniques.

Intended Use

Bioptimal Angiographic Kits are designed for use in the cath lab during PTCA procedures or Angiography to monitor patient's arterial pressure, administer intravenous solution and contrast media to the patient.

Angiographic Kit Components

- Multi-Gang Manifold
- IV Administration Set
- Contrast Media Administration Set
- Control Syringes
- High Pressure tubing
- Pressure Tubing
- IV Tubing

Angiographic Kits		
TA2000	Prime Angiographic Kit II	2-GMF, RH(OFF)-500PSI
TA2001	3-GMF, RH(OFF)-200PSI	DPT-AMG#3039
TA3000	3-GMF, RH(OFF)-550PSI	CGH 3-Gang Manifold Kit
TA3001	3-GMF, RH(ON)-550PSI	CGH 3-Gang Manifold Kit W/4 Syringes
CMC-332	2-GMF, RH(ON)-500PSI	Angiographic Kit
BMB Angio Kit	3-GMF, RH(OFF)-500PSI	
Prime Angiographic Kit	3-GMF, RH(ON)-500PSI	