



---

# **ASAHI CHIKAI™ FAMILY - PRODUCT DESCRIPTION -**



## ASAHI CHIKAI™

- ASAHI CHIKAI™ is designed as workhorse guide wire for neuro interventional procedures. CHIKAI family are as follows;
  - ✓ CHIKAI™
  - ✓ CHIKAI™ 10
  - ✓ CHIKAI™ 008
- CHIKAI and CHIKAI 10 are designed for aneurysm treatment and CHIKAI 008 is designed to be used with flow guided microcatheters.


ASAHI  
**CHIKAI**  
Neurovascular GuideWire

ASAHI  
**CHIKAI 10**  
Neurovascular GuideWire

ASAHI  
**CHIKAI 008**  
Neurovascular GuideWire




# Initial Benchmarks of ASAHI CHIKAI™ Guide Wire





# ASAHI CHIKAI™ Structure




- **Usable Length : 200cm, 300cm**
- **Core Wire Material: Stainless Steel**
- **Distal/Proximal OD : 0.014" (0.36mm)**
- **Hydrophilic Coating : 170cm from Distal**
- **PTFE Coating : 130cm from proximal (300cm only)**
- **Spring Coil Length : 30cm**
- **Radiopaque Length : 5cm**



# ASAHI CHIKAI™ 10 Structure

ASAHI  
CHIKAI 香 AI 10  
Neurovascular GuideWire



- Usable Length : 200cm, 300cm
- Core Wire Material: Stainless Steel
- Distal/Proximal OD : 0.010" (0.26mm)
- Hydrophilic Coating : 170cm from Distal
- PTFE Coating : 130cm from proximal (300cm only)
- Spring Coil Length : 9.5cm
- Radiopaque Length : 3cm



## ASAHI CHIKAI™ Key Features

---

- Torqueability
  - ✓ Able to negotiate with tortious vessels
  - ✓ Able to select branches
- Deliverability and lubricity
  - ✓ Keep tip shape and deliver microcatheter precisely
- Conformability
  - ✓ Not stretch vessels
- Tip shape memory/Durability
  - ✓ Good shape retention
- Variety of size
  - ✓ 300cm is available for microcatheter exchange
  - ✓ CHIKAI™ 200cm is connectable to ASAHI Extension Wire