

# Figulla® Flex II The Third Generation ASD and PFO Occluders

# Compatibility Chart

## Figulla® Flex II PFO Occluder (Single LA Disc)

Article No.	Product Name	Disc Ø (mm)	Sheath Size (F)
18PFO25S	Figulla® Flex II PFO	23/25	9

## Figulla® Flex II PFO Occluder

Article No.	Product Name	Disc Ø (mm)	Sheath Size (F)
19PFO18D	Figulla® Flex II PFO	16/18	7
19PFO25D	Figulla® Flex II PFO	23/25	9
19PFO30D	Figulla® Flex II PFO	27/30	9
19PFO35D	Figulla® Flex II PFO	31/35	11


## Figulla® Flex II ASD Occluder

Article No.	Product Name	Waist Ø (mm)	Sheath Size (F)
29ASD04*	Figulla® Flex II ASD	Size 4	7
29ASD05*	Figulla® Flex II ASD	Size 5	7
29ASD06	Figulla® Flex II ASD	Size 6	7
29ASD07	Figulla® Flex II ASD	Size 7.5	7
29ASD09	Figulla® Flex II ASD	Size 9	7
29ASD10	Figulla® Flex II ASD	Size 10.5	7
29ASD12	Figulla® Flex II ASD	Size 12	9
29ASD13*	Figulla® Flex II ASD	Size 13.5	9
29ASD15	Figulla® Flex II ASD	Size 15	9
29ASD16*	Figulla® Flex II ASD	Size 16.5	9
29ASD18	Figulla® Flex II ASD	Size 18	9
29ASD21	Figulla® Flex II ASD	Size 21	11
29ASD24	Figulla® Flex II ASD	Size 24	11
29ASD27	Figulla® Flex II ASD	Size 27	12
29ASD30	Figulla® Flex II ASD	Size 30	12
29ASD33	Figulla® Flex II ASD	Size 33	12
29ASD36	Figulla® Flex II ASD	Size 36	12
29ASD39	Figulla® Flex II ASD	Size 39	12
29ASD40	Figulla® Flex II ASD	Size 40	12

## Figulla® Flex II Uniform Occluder

Article No.	Product Name	Disc Ø (mm)	Sheath Size (F)
16UNI17	Figulla® Flex II UNI	17/17	7
16UNI24	Figulla® Flex II UNI	24/24	9
16UNI28	Figulla® Flex II UNI	28.5/28.5	9
16UNI33	Figulla® Flex II UNI	33/33	11
16UNI40*	Figulla® Flex II UNI	40/40	12

\* Subject to availability




# Figulla® Flex II The Third Generation ASD and PFO Occluders


## Occlutech® Delivery Set

Product Range - Dimensions and sizes.

Art No.	Size (F)	Length Sheath (mm)	Length Loader (mm)	Dilator Tip (mm)	Hemo-valve (mm)	Distal Angle
51DS007	7	800	120	25	70	45°
51DS008	8	800	160	25	70	45°
51DS009	9	800	160	25	70	45°
51DS010	10	800	160	25	70	45°
51DS011	11	800	200	25	70	45°
51DS012	12	800	200	25	70	45°
51DS014	14	800	200	25	70	45°



The Flex II braid allows on average **20 % less diameter in the delivery sheath**



## Perfecting Performance

### Customer Service International

Occlutech, Rönnowsgatan 8, S-252 25 Helsingborg, Sweden  
Tel +46 (0)42 33 65 25, Fax +46 (0)42 311 09 70  
order@occlutech.com, [www.occlutech.com](http://www.occlutech.com)



© Occlutech 2012. All rights reserved, P07F06.017.01

For marketing purposes only. Not subject to update. For updated information please see instructions for use. (P07F05.013\_02)