


# ASAHI Gaia

P T C A G U I D E W I R E


# Ordering Information

## ■ Structure



## ■ Composite Core

Multiple wire components to enhance durability and torque transmission.



Employs the ASAHI brand rope coil,  
 Which provides torqueability, flexibility and resiliency.

Product	Catalog No.	Diameter	Coating Type	Usable Length	Spring Coil Length	Radiopaque Length	Tip Shape
ASAHI Gaia First	AHW14R007P	0.36 / 0.26mm (0.014 / 0.010 inch)	SLIP-COAT®	190cm	15cm	15cm	1mm pre-shape
	AHW14R307P			300cm			
ASAHI Gaia Second	AHW14R008P	0.36 / 0.28mm (0.014 / 0.011 inch)	SLIP-COAT®	190cm	15cm	15cm	1mm pre-shape
	AHW14R308P			300cm			
ASAHI Gaia Third	AHW14R011P	0.36 / 0.30mm (0.014 / 0.012 inch)	SLIP-COAT®	190cm	15cm	15cm	1mm pre-shape
	AHW14R311P			300cm			

Coated with SLIP-COAT® coating.  
 SLIP-COAT® is a trademark of Angiotech BioCoatings Corp. in the U.S. and /or other countries.

ASAHI Gaia is a trademark or registered trademark of ASAHI INTECC CO., LTD. in Japan and other countries.



Manufactured by:



**ASAHI INTECC CO.,LTD.**

3-100 Akatsuki-cho, Seto, Aichi 489-0071 Japan  
 TEL:+81-561-48-5551 FAX:+81-561-48-5552  
<http://www.asahi-intecc.com>

Distributed by: