

Hemostasis Valve Kit (Push-Click)

Patent design of profile for ease of use
 One step valve control: Click-Open/Click-Close
 Valve can be adjusted between Open/Close
 400Kpa leak-proof seal from proprietary valve
 Fully open for 9 French
 2 types of extension for option
 Torquer and Introducer fits to 0.014" – 0.021" guidewire
 Various kit configuration available



CE₀₁₂₃ **ARTG** **SFDA**

Ref No.	Bypass Extension	Guide Wire Torque Device	Guide Wire Insertion Tools
DPY01	N	N	N
DPY01B		Torquer (Plastic Buckle)	Y
DPY01R		Torquer (Metal Buckle)	
DPY02	Fixed Extension Tube & Stopcock	N	N
DPY02B		Torquer (Plastic Buckle)	Y
DPY02R		Torquer (Metal Buckle)	
DPY03	Removable Extension Tube & Stopcock	N	N
DPY03B		Torquer (Plastic Buckle)	Y
DPY03R		Torquer (Metal Buckle)	

International Marketing

Tel: +61 2 800 617 46

Fax: +61 2 808 864 77

Email: marketing@demaxmedical.com