

# Thrombuster II

## Features

- 1**
  - "Big" and "round" suction lumen
  - Minimum flow resistance enhances the suction performance
- 2**
  - Shortened wire lumen length enhances flexibility of the distal shaft
- 3**
  - Hydrophilic coating (30 cm from the tip) assures smooth insertion into the blood vessel
- 4**
  - Removable core wire enhances pushability and anti-kinking feature



## Product specifications

Catalogue Number	Minimum Inner Diameter of Compatible GC	Compatible GW	Shaft Diameter				Guidewire Lumen Part		Effective Length
			Outer		Inner		Maximum Outer Diameter	Length	
			Distal	Proximal	Distal	Proximal			
T2R6SA4W	0.070inch (1.78mm)	0.014inch (0.36mm)	3.9Fr (1.30mm)	3.9Fr (1.30mm)	3.0Fr (1.00mm)	3.3Fr (1.10mm)	5.1Fr (1.70mm)	10mm	1400mm
T2R7SA4W	0.080inch (2.03mm)	0.014inch (0.36mm)	4.5Fr (1.50mm)	4.6Fr (1.53mm)	3.6Fr (1.20mm)	3.9Fr (1.30mm)	5.7Fr (1.90mm)	10mm	1400mm




For single use only  
Do not reuse or resterilize

STERILE EO


Sterilized by ethylene oxide gas



Read the accompanying instructions for use carefully



Side View of Distal Tip



Manufacturer

**KANEKA CORPORATION**

3-2-4, Nakanoshima, Kita-ku, Osaka, Japan  
TEL. +81-6-6226-5285 FAX. +81-6-6226-4506

EC-Representative

**KANEKA PHARMA EUROPE N.V.**

Triomflaan 173, B-1160 Brussels, Belgium  
TEL. +32-(0)-2-663-0340 FAX. +32-(0)-2-663-0361

