

UniQual


Проводник для ангиографии с нитиноловым стержнем

UniQual *Успех манипуляции достигается его стабильными и надёжными характеристиками*

Гибкий кончик проводника

Гибкий кончик проводника минимизирует травму сосуда и способствует упрощению хода манипуляции


Тест нагрузки кончика



Отличная проталкиваемость

Высокая жёсткость центрального стержня для лучшей управляемости и проталкиваемости в процессе введения проводника

Тест жёсткости



Выраженная рентгеноконтрастность

Выраженная рентгеноконтрастность проводника при флюороскопии упрощает манипуляции.


Тест рентгеноконтрастности



Информация для заказа

Номер по каталогу	Наружный диаметр	Форма кончика	Длина (см)
US0035B1500	0.035	Прямой	150
US0035B2600			260
UA0035B0500	0.035	Изогнутый под углом	50
UA0035B0800			80
UA0035B1500			150
UA0035B1800			180
UA0035B2000			200
UA0035B2200			220
UA0035B2600			260
UA0035B3000			300
UJ1535B1500	0.035	1.5J-Tip	150
UJ1535B1800			180
UJ1535B2000			200
UJ1535B2200			220
UJ1535B2600			260
UJ1535B3000			300

Номер по каталогу	Наружный диаметр	Форма кончика	Длина (см)
UJ3035B1500	0.035	3J-Tip	150
UJ3035B2000			200
UJ3035B2600			260
UJ3035B3000			300
UJ6035B1500	0.035	6J-Tip	150
UJ6035B2000			200
UJ6035B2600			260
UJ6035B3000			300



Изготовитель:



ASAHI INTECC CO.,LTD.

3-100 Akatsuki-cho, Seto, Aichi 489-0071 Japan
 TEL: +81-561-48-5551 FAX: +81-561-48-5552
<http://www.asahi-intecc.com>

Дистрибьютор:



ASAHI INTECC CO., LTD. EUROPE

Strawinskylaan 967
 (WTC Tower D-9th floor)
 1077 XX Amsterdam,
 The Netherlands
 TEL: +31-(0)20-794-0640
 FAX: +31-(0)20-794-0641

ASAHI INTECC USA, INC.

2500 Red Hill Avenue, Suite 210
 Santa Ana, California,
 92705, USA
 TEL: +1-949-756-8901
 FAX: +1-949-756-8165

ASAHI INTECC (THAILAND) CO., LTD.

158/1 Moo5,
 Bangkadi Industrial Park,
 Tiwanon Road, Tambol Bangkadi,
 Amphur Muang,
 Pathumthani 12000,
 Thailand
 TEL: +66-2-501-1302
 FAX: +66-2-501-1305

ASAHI INTECC CO., LTD.
 SINGAPORE BRANCH

3 Bishan Place
 #06-01 Bishan CPE Building
 Singapore 579838
 TEL: +65-6250-4248
 FAX: +65-6250-4247

ASAHI INTECC CO., LTD.

BEIJING REPRESENTATIVE OFFICE
 Room 18, 10th Floor,
 Twin Towers (East),
 B12 Jianguomenwai Ave.,
 Beijing, China 100022
 TEL: +86-10-5123-5100
 FAX: +86-10-5123-5177

ASAHI INTECC CO., LTD.
 HONG KONG BRANCH

Unit 501-2, 5/F, Hon Kwok Jordan Centre
 No 7 Hillwood Road,
 Tsim Sha Tsui, Kowloon
 Hong Kong
 TEL: +852-2316-2510
 FAX: +852-2316-2485